



CoreLink.
The Source for Spine™

INSIGHT | PERFORMANCE | VALUE

FLXfit[®] 15


ARTICULATING EXPANDABLE POSTERIOR LUMBAR SPACER

EXPANDABLE

Up to 15° of controlled
and continuous
lordotic expansion


ARTICULATING

Integrated articulation
for optimal anterior
midline placement



THE FLXfit® 15 ADVANTAGE



FLXfit15 is an expandable and articulating posterior lumbar interbody that allows for controlled, in-situ expansion.



SYSTEM FEATURES

- Up to 4mm of anterior height expansion
- Up to 15° of controlled and continuous lordotic expansion
- Open device architecture enables graft material packing
- Bullet-nosed to optimize insertion
- Inferior and superior surfaces include teeth intended to reduce cage migration
- Single instrument for insertion and expansion of cage

FLXfit15 FOOTPRINT DIMENSIONS



CAGE AND CARTRIDGE PART NUMBERS

CATALOG #	DESCRIPTION
AT3208	32MM x 9MM x 8MM
AT3209	32MM x 9MM x 9MM
AT3210	32MM x 9MM x 10MM
AT3211	32MM x 9MM x 11MM
AT3212	32MM x 9MM x 12MM
AT3213	32MM x 9MM x 13MM

CATALOG #	DESCRIPTION
AT4008	40MM x 9MM x 8MM
AT4009	40MM x 9MM x 9MM
AT4010	40MM x 9MM x 10MM
AT4011	40MM x 9MM x 11MM
AT4012	40MM x 9MM x 12MM
AT4013	40MM x 9MM x 13MM
SD-15-FU	FLXfit15 CARTRIDGE (SINGLE USE)



Manufactured in the
United States of America.

St. Louis, MO
CoreLinkSurgical.com

P: (888) 349-7808
F: (888) 629-5715

CL-MKTG-0073 Rev. 2